HAPILIN[®]

ULTRASONIC SKIN SCRUBBER

www.hapilin.com

HAPILIN[®]

ULTRASONIC SKIN SCRUBBER

www.hapilin.com

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Letter of Guarantee:

This appliance has a CE and ROHS marking that complies to the health, safety, and environmental protection standards of the European Economic Area. It also complies to the American Quality and Safety Standards.

It comes with a 14-day money-back guarantee which starts from the day of invoice on the receipt. It can extend to 30 days in the case of a manufacturing defect. It also comes with a 2-year repair warranty.

In the event your appliance stops working within the first 2 years of purchase, please contact us by sending an email to warranty@hapilin.com explaining the problem clearly. We will then be happy to advise you on what to do to ensure your satisfaction.

The warranty applies to manufacturing defects only and does not apply to the following cases:

- Negligence and product misuse.
- Breaking or cracking the appliance.
- Exposing the appliance to water.
- Damage caused by an attempt to repair the appliance by an unauthorized party.

Important: The money-back guarantee is not applied in the case of personal use of the product.

The appliance has been inspected prior to shipment, so any damage or crack occurring to the outer casing is not covered by this warranty.

As a company, we will be responsible for any costs incurred for repair/replacement of your appliance in the event of a manufacturing fault, in addition to the return shipping costs. However, you will be responsible for shipping costs for the return of the original appliance, if the fault is not a manufacturing one, which must be sent back using a trackable/traceable delivery service, returned in its original box with the original purchase receipt.

Please keep this warranty document with your appliance so that it is handy in the event a problem arises.

Cleaning:

Rinse the scrubber probe and wipe it with a dry soft cloth after using.

Warning: Do not immerse the appliance into water. Only clean it with a dry cloth from the outside.

Safety Matters:

- Pregnant women, diabetic patients, and people with heart diseases are forbidden to use this appliance.
- Make sure to consult a doctor before using this appliance if you have any filler injections.
- People with skin inflammation, scars, or electromagnetic sensitivity is forbidden to use this appliance.
- Always inspect the appliance before use. Do not use if damaged.
- Don't wear fiber clothes before using.
- Take off watches, glasses and other metal objects before using.
- Make sure to use with dry hands.
- Avoid using the scrubber probe close to eyeball.
- Don't operate in the same skin area for a long time.
- Always clean the appliance after each use.
- Avoid using acetone, alcohol, benzene, or any other abrasive cleansers that could dull, bleach, or harm the plastic.
- This appliance should not be used, cleaned, or maintained by children.
- Do not use on children below 14 years without adult supervision.
- Keep it away from high temperature.
- Please do not reassemble, modify, or deconstruct the appliance in any way.

EMS/H: Turn on the device into SONIC & EMS/H by pressing 4 times on the ON/OFF button.

Manual Shut-Off: Switch back to EMS/H mode, then press the ON/OFF button to turn it off.

Note: The device shuts off automatically after 5 minutes of use.

Usage Instructions:

1. Make sure to clean your face before using.

2. Apply some skin lotion, cream, gel, or serum to avoid using it on dry skin.

3. Turn on the (SONIC & ION+) mode for deep-cleansing, by gently gliding the metal probe edge along the skin's surface while facing downwards.

4. Turn on the (SONIC & ION-) and (SONIC & EMS) for improving nutrient absorption and skin tightening, by gently gliding the metal probe surface on the skin.

5. Watch as the gunk from the pores easily builds up onto the edge of the probe.

6. Wash your face with water after treatment.

7. It is recommended to circulate your face with ice cubes to close face pores.

Note: It is recommended to use this appliance 2-3 times/week.

Product Features:

Cordless and Rechargeable:

It can be charged by placing it on the charging base, then connect it to a suitable electrical outlet.

One-Touch Operation:

The appliance has only one button for easy using and fast shifting between modes.

Ergonomic Design:

Ergonomically curved to fit comfortably in your hand. It gives you maximum comfort, efficiency, and safety to use.

LED Light Indicators:

- **Automatic Shut-Off:** The appliance beeps 2 times, then LED light shuts off.
- **Low Battery:** Beeps 3 times, LED light keeps flashing in red.
- **In Charging:** LED light keeps flashing in red.
- **Fully Charged:** LED light turns green.

Operating Instructions:

Ultrasonic Exfoliation (SONIC & ION+):

Turn on the device into SONIC & ION+ by pressing once on the ON/OFF button.

Anion Lead-In (SONIC & ION-):

Turn on the device into SONIC & ION- by pressing twice on the ON/OFF button.

EMS/L:

Turn on the device into SONIC & EMS/L by pressing 3 times on the ON/OFF button.

Specifications:

Product Name: Hapilin® Ultrasonic Skin Scrubber

Model No. HA-1610

Rated Voltage: 5V

Rated Power: 2.5W

Rated Current: ≤ 500 mA

Power Supply: Inductive Charging Base

Ultrasonic Frequency: 27 ± 5 kHz

Charging Time: ≤ 4 hrs

Operating Time: ≥ 50 mins

Timer: 5 mins

Product Functions:

Ultrasonic Exfoliation (SONIC & ION+):

It is a skin cleansing procedure that creates ultrasonic waves and high-frequency vibration that help dead skin cell removal. It also boosts collagen production, promotes cellular renewal and increases blood circulation. In addition, the blue photon light gives a double cleaning effect, cleans pores, prevents bacterial multiplication and balances sebum secretion.

Anion Lead-In (SONIC & ION-):

The anion attracts human micro-electricity, as the metal scrubber probe effectively leads the nutrition into the deep skin layer. In addition, the red photon light with penetrability, improves the absorption of nutrients, leaving you nourished with brightened skin.

Electrical Muscle Stimulation EMS (SONIC & EMS):

The circulating current formed by the metal scrubber probe restores skin elasticity that helps skin tightening, firming and smoothing.

Hapilin® Ultrasonic Skin Scrubber

Hapilin® Ultrasonic Skin Scrubber is a multifunctional beauty device that features a stainless steel probe that provides ultrasonic oscillation technology which gently exfoliates the skin and removes oil, dirt, and debris. It also features +/- ionic reaction and EMS technologies. It is a 3-step routine for regaining natural beauty, as it:

- ▶ Provides unique therapeutic massage and high cellular metabolism.
- ▶ Cleanses clogged pores.
- ▶ Helps in deep cleansing and skin rejuvenation.
- ▶ Promotes blood circulation through micro-patting.
- ▶ Extracts blackheads and whiteheads.
- ▶ Exfoliates dry and flaky skin making your skin smooth and delicate.
- ▶ Helps in skin tightening and lifting.
- ▶ Removes acne scars and eliminates the appearance of skin pores.
- ▶ Better skin absorption of skincare products.

- ① Protective Cover
- ② Stainless Steel Scrubber Probe
- ③ Photon Light (Red-Blue)
- ④ Ultrasonic Indicator
- ⑤ Positive Ion Indicator
- ⑥ Negative Ion (Anion) Indicator
- ⑦ Low EMS Indicator (EMS/L)
- ⑧ High EMS Indicator (EMS/H)
- ⑨ ON/OFF Button
- ⑩ Inductive Charging Base
- ⑪ Universal Voltage Charger